

For Research Use Only

TFDP2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84408-2-PBS



Basic Information

| | | |
|--|--|---|
| Catalog Number: 84408-2-PBS | GenBank Accession Number: BC021113 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 7029 | CloneNo.: 241457A11 |
| Source: Rabbit | UNIPROT ID: Q14188 | |
| Isotype: IgG | Full Name: transcription factor Dp-2 (E2F dimerization partner 2) | |
| Immunogen Catalog Number: AG2054 | Calculated MW: 446 aa, 49 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01288-1, TFDP2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84408-1-PBS. Detection antibody: 84408-2-PBS. Standard: Ag2054. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP01288-3, TFDP2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84408-3-PBS. Detection antibody: 84408-2-PBS. Standard: Ag2054. Range: 0.313-40 ng/mL