

For Research Use Only

MLL4 Recombinant antibody

Catalog Number: 84410-1-RR



Basic Information

Catalog Number:

84410-1-RR

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG35385

GenBank Accession Number:

BC009337

GeneID (NCBI):

9757

UNIPROT ID:

Q9UMN6

Full Name:

myeloid/Lymphoid or mixed-Lineage
leukemia 4

Observed MW:

294 kDa

Purification Method:

Protein A purification

CloneNo.:

241819A4

Recommended Dilutions:

WB 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, rat

Positive Controls:

WB : PC-12 cells, HEK-293T cells, A549 cells

Background Information

MLL4, also known as KMT2B, which is an essential subunit of the H3K4 methylation complexes, regulates cell type-specific transcriptional programs through enhancer activation. Taspase 1 mediates cleavage of the MLL protein, generating cleavage products.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 84410-1-RR (MLL4 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



Biolayer interferometry (BLI) kinetic assays of 84410-1-RR against Human MLL4 were performed. The affinity constant is 22.6 pM.