

For Research Use Only

IL1A Recombinant antibody, PBS Only (Detector)

Catalog Number: 84411-1-PBS



Basic Information

Catalog Number: 84411-1-PBS	GenBank Accession Number: BC013142	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 3552	CloneNo.: 241752B8
Source: Rabbit	UNIPROT ID: P01583	
Isotype: IgG	Full Name: interleukin 1, alpha	
Immunogen Catalog Number: HZ-1320	Calculated MW: 271 aa, 31 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01289-1, IL-1 alpha Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84411-2-PBS. Detection antibody: 84411-1-PBS. Standard: HZ-1320. Range: 0.039-10 ng/mL.