

For Research Use Only

# CCL20/MIP-3 alpha Recombinant antibody, PBS Only (Detector)

Catalog Number: 84413-2-PBS



## Basic Information

<b>Catalog Number:</b> 84413-2-PBS	<b>GenBank Accession Number:</b> NM_004591.3	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 6364	<b>CloneNo.:</b> 241694B2
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> P78556-1	
<b>Isotype:</b> IgG	<b>Full Name:</b> chemokine (C-C motif) ligand 20	
	<b>Calculated MW:</b> 11kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01301-2, CCL20/MIP-3 alpha Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84413-1-PBS. Detection antibody: 84413-2-PBS. Standard: Eg2398. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP01301-3, CCL20/MIP-3 alpha Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84413-3-PBS. Detection antibody: 84413-2-PBS. Standard: Eg2398. Range: 0.781-100 ng/mL