

For Research Use Only

IL-4 R alpha/CD124 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84414-1-PBS



Basic Information

Catalog Number:	84414-1-PBS	GenBank Accession Number:	NM_000418.4	Purification Method:	Protein A purification
Size:	1 mg/ml	GeneID (NCBI):	3566	CloneNo.:	241824D2
Source:	Rabbit	UNIPROT ID:	P24394-1		
Isotype:	IgG	Full Name:	interleukin 4 receptor		
		Calculated MW:	90kDa		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01282-3, IL-4 R alpha/CD124 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84422-5-PBS. Detection antibody: 84422-2-PBS. Standard: Eg1684. Range: 0.156-20 ng/mL