

For Research Use Only

IL-4R/CD124 Recombinant antibody, PBS Only

Catalog Number: 84414-6-PBS



Basic Information

Catalog Number:

84414-6-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_000418.4

GeneID (NCBI):

3566

UNIPROT ID:

P24394-1

Full Name:

interleukin 4 receptor

Calculated MW:

90 kDa

Purification Method:

Protein A purification

CloneNo.:

241824A11

Applications

Tested Applications:

IF/ICC, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed Daudi cells using IL-4R/CD124 antibody (84414-6-RR, Clone: 241824A11) at dilution of 1:250 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 84414-6-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 84414-6-RR against Human IL-4R/CD124 were performed. The affinity constant is below 1 pM.