

For Research Use Only

ZNHIT2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84418-4-PBS



Basic Information

Catalog Number: 84418-4-PBS	GenBank Accession Number: BC052240	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 741	CloneNo.: 241868C11
Source: Rabbit	UNIPROT ID: Q9UHR6	
Isotype: IgG	Full Name: zinc finger, HIT type 2	
Immunogen Catalog Number: AG10351	Calculated MW: 403 aa, 43 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01304-1, ZNHIT2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84418-3-PBS. Detection antibody: 84418-4-PBS. Standard: Ag10351. Range: 0.781-100 ng/mL.