

For Research Use Only

BCMA/TNFRSF17 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84419-2-PBS



Basic Information

Catalog Number:

84419-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC058291

GeneID (NCBI):

608

UNIPROT ID:

Q02223

Full Name:

tumor necrosis factor receptor
superfamily, member 17

Calculated MW:

184 aa, 20 kDa

Purification Method:

Protein A purification

CloneNo.:

241731A11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01305-1, TNFRSF17 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84419-4-PBS. Detection antibody: 84419-2-PBS. Standard: Eg1812. Range: 0.313-20 ng/mL.