

For Research Use Only

TNFRSF17 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84419-2-PBS



Basic Information

| | | |
|---------------------------------------|---|---|
| Catalog Number: 84419-2-PBS | GenBank Accession Number: BC058291 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 608 | CloneNo.: 241731A11 |
| Source: Rabbit | UNIPROT ID: Q02223 | |
| Isotype: IgG | Full Name: tumor necrosis factor receptor superfamily, member 17 | |
| | Calculated MW: 184 aa, 20 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01305-1, TNFRSF17 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84419-4-PBS. Detection antibody: 84419-2-PBS. Standard: Eg1812. Range: 0.313-20 ng/mL.