

For Research Use Only

# TNFRSF17 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84419-4-PBS



## Basic Information

Catalog Number:

84419-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC058291

GeneID (NCBI):

608

UNIPROT ID:

Q02223

Full Name:

tumor necrosis factor receptor superfamily, member 17

Calculated MW:

184 aa, 20 kDa

Purification Method:

Protein A purification

CloneNo.:

241731E7

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01305-1, TNFRSF17 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84419-4-PBS. Detection antibody: 84419-2-PBS. Standard: Eg1812. Range: 0.313-20 ng/mL.