

For Research Use Only

Rat IL-4 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84422-2-PBS



Basic Information

Catalog Number:

84422-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_201270.1

GeneID (NCBI):

287287

UNIPROT ID:

P20096

Full Name:

interleukin 4

Calculated MW:

16kDa

Purification Method:

Protein A purification

CloneNo.:

241249D5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Cytometric bead array standard curve of MP01282-2, IL-4 R alpha/CD124 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84422-3-PBS. Detection antibody: 84422-2-PBS. Standard: Eg1113. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP01282-3, IL-4 R alpha/CD124 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84422-5-PBS. Detection antibody: 84422-2-PBS. Standard: Eg1113. Range: 1.25-80 ng/mL



Cytometric bead array standard curve of MP01293-2, RAT IL-4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84422-3-PBS. Detection antibody: 84422-2-PBS. Standard: Eg1113. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP01293-3, RAT IL-4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84422-5-PBS. Detection antibody: 84422-2-PBS. Standard: Eg1113. Range: 1.25-80 ng/mL