

For Research Use Only

Rat IL-4 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84422-4-PBS



Basic Information

Catalog Number:	84422-4-PBS	GenBank Accession Number:	NM_201270.1	Purification Method:	Protein A purification
Size:	1 mg/ml	GeneID (NCBI):	287287	CloneNo.:	241249G7
Source:	Rabbit	UNIPROT ID:	P20096		
Isotype:	IgG	Full Name:	interleukin 4		
		Calculated MW:	16kDa		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
rat

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01282-1, IL-4 R alpha/CD124 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84422-4-PBS. Detection antibody: 84422-1-PBS. Standard: Eg1113. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP01293-1, RAT IL-4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84422-4-PBS. Detection antibody: 84422-1-PBS. Standard: Eg1113. Range: 0.625-80 ng/mL