

For Research Use Only

# HSPA13 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number:84426-1-PBS



## Basic Information

<b>Catalog Number:</b> 84426-1-PBS	<b>GenBank Accession Number:</b> BC036370	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 6782	<b>CloneNo.:</b> 241572B3
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> P48723	
<b>Isotype:</b> IgG	<b>Full Name:</b> heat shock protein 70kDa family, member 13	
<b>Immunogen Catalog Number:</b> AG3349	<b>Calculated MW:</b> 471 aa, 52 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01302-1, HSPA13 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84426-1-PBS. Detection antibody: 84426-2-PBS. Standard: Ag3349. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP01302-3, HSPA13 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84426-3-PBS. Detection antibody: 84426-1-PBS. Standard: Ag3349. Range: 0.781-100 ng/mL