

For Research Use Only

Mouse TNFSF11 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84431-1-PBS



Basic Information

Catalog Number: 84431-1-PBS	GenBank Accession Number: NM_011613.3	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 21943	CloneNo.: 241791A7
Source: Rabbit	UNIPROT ID: O35235	
Isotype: IgG	Full Name: tumor necrosis factor (ligand) superfamily, member 11	
	Calculated MW: 35KD	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01276-1, MOUSE TNFSF11/RANKL Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84431-1-PBS. Detection antibody: 84431-4-PBS. Standard: Eg1956. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP01276-3, MOUSE TNFSF11/RANKL Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84431-1-PBS. Detection antibody: 84431-2-PBS. Standard: Eg1956. Range: 0.313-40 ng/mL