

For Research Use Only

Mouse TNFSF11 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84431-3-PBS



Basic Information

Catalog Number:

84431-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_011613.3

GeneID (NCBI):

21943

UNIPROT ID:

O35235

Full Name:

tumor necrosis factor (ligand)
superfamily, member 11

Calculated MW:

35KD

Purification Method:

Protein A purification

CloneNo.:

241791B11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01276-2, MOUSE TNFSF11/RANKL Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84431-3-PBS. Detection antibody: 84431-4-PBS. Standard: Eg1956. Range: 0.313-40 ng/mL.