

For Research Use Only

ATG9A Recombinant antibody, PBS Only (Detector)

Catalog Number: 84433-2-PBS



Basic Information

Catalog Number: 84433-2-PBS	GenBank Accession Number: BC065534	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 79065	CloneNo.: 241652F9
Source: Rabbit	UNIPROT ID: Q7Z3C6	
Isotype: IgG	Full Name: ATG9 autophagy related 9 homolog A (S. cerevisiae)	
Immunogen Catalog Number: AG24212	Calculated MW: 837 aa, 94 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01278-1, ATG9A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84433-1-PBS. Detection antibody: 84433-2-PBS. Standard: Ag24212. Range: 0.313-40 ng/mL