

For Research Use Only

Mouse 5'-Nucleotidase/CD73 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84434-5-PBS



Basic Information

Catalog Number: 84434-5-PBS	GenBank Accession Number: NM_011851.4	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 23959	CloneNo.: 241541G11
Source: Rabbit	UNIPROT ID: Q61503	
Isotype: IgG	Full Name: 5' nucleotidase, ecto	
	Calculated MW: 64kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01279-2, MOUSE CD73 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84434-5-PBS. Detection antibody: 84434-2-PBS. Standard: Eg1389. Range: 0.625-40 ng/mL.



Cytometric bead array standard curve of MP01279-3, MOUSE CD73 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84434-5-PBS. Detection antibody: 84434-1-PBS. Standard: Eg1389. Range: 0.625-40 ng/mL.