

For Research Use Only

TANC2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84435-5-PBS



Basic Information

Catalog Number:

84435-5-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG35427

GenBank Accession Number:

NM_001394998

GeneID (NCBI):

26115

UNIPROT ID:

Q9HCD6

Full Name:

tetratricopeptide repeat, ankyrin
repeat and coiled-coil containing 2

Calculated MW:

219kDa

Purification Method:

Protein A purification

CloneNo.:

241631E7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01281-1, TANC2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84435-2-PBS. Detection antibody: 84435-5-PBS. Standard: Ag35427. Range: 0.781-50 ng/mL



Cytometric bead array standard curve of MP01281-2, TANC2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84435-1-PBS. Detection antibody: 84435-5-PBS. Standard: Ag35427. Range: 0.781-100 ng/mL