

For Research Use Only

Mouse Myeloperoxidase Recombinant antibody, PBS Only (Detector)

Catalog Number: 84436-1-PBS



Basic Information

| | | |
|--------------------------------|---------------------------------------|--|
| Catalog Number: 84436-1-PBS | GenBank Accession Number: AAR99349 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 17523 | CloneNo.: 241839A11 |
| Source: Rabbit | UNIPROT ID: P11247 | |
| Isotype: IgG | Full Name: myeloperoxidase | |
| | Calculated MW: 81KD | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01284-1, MOUSE MPO Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84436-2-PBS. Detection antibody: 84436-1-PBS. Standard: Eg0702. Range: 1.25-80 ng/mL



Cytometric bead array standard curve of MP01284-2, MOUSE MPO Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84436-3-PBS. Detection antibody: 84436-1-PBS. Standard: Eg0702. Range: 1.25-80 ng/mL