

For Research Use Only

Mouse Myeloperoxidase Recombinant antibody, PBS Only (Detector)

Catalog Number: 84436-1-PBS



Basic Information

Catalog Number:

84436-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

AAR99349

GeneID (NCBI):

17523

UNIPROT ID:

P11247

Full Name:

myeloperoxidase

Calculated MW:

81KD

Purification Method:

Protein A purification

CloneNo.:

241839A11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01284-1, MOUSE MPO Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84436-2-PBS. Detection antibody: 84436-1-PBS. Standard: Eg0702. Range: 1.25-80 ng/mL



Cytometric bead array standard curve of MP01284-2, MOUSE MPO Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84436-3-PBS. Detection antibody: 84436-1-PBS. Standard: Eg0702. Range: 1.25-80 ng/mL