

For Research Use Only

Mouse Myeloperoxidase Recombinant antibody, PBS Only (Capture)

Catalog Number: 84436-2-PBS



Basic Information

Catalog Number: 84436-2-PBS	GenBank Accession Number: AAR99349	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 17523	CloneNo.: 241839E2
Source: Rabbit	UNIPROT ID: P11247	
Isotype: IgG	Full Name: myeloperoxidase	
	Calculated MW: 81KD	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01284-1, MOUSE MPO Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84436-2-PBS. Detection antibody: 84436-1-PBS. Standard: Eg0702. Range: 1.25-80 ng/mL



Cytometric bead array standard curve of MP01284-3, MOUSE MPO Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84436-2-PBS. Detection antibody: 84436-3-PBS. Standard: Eg0702. Range: 0.625-80 ng/mL