

For Research Use Only

Mouse Nectin-2/CD112 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84439-1-PBS



Basic Information

Catalog Number: 84439-1-PBS	GenBank Accession Number: NM_008990.3	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 19294	CloneNo.: 241798A2
Source: Rabbit	UNIPROT ID: P32507-1	
Isotype: IgG	Full Name: poliovirus receptor-related 2	
	Calculated MW: 57kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01292-1, MOUSE Nectin-2/CD112 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84439-3-PBS. Detection antibody: 84439-1-PBS. Standard: Eg2019. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP01292-2, Mouse Nectin-2/ CD112 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84439-3-PBS. Detection antibody: 84439-1-PBS. Standard: Eg2019. Range: 0.313-20 ng/mL