

For Research Use Only

ADAP1 Recombinant antibody

Catalog Number: 84443-1-RR



Basic Information

Catalog Number:

84443-1-RR

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4898

GenBank Accession Number:

BC033747

GeneID (NCBI):

11033

UNIPROT ID:

O75689

Full Name:

ArfGAP with dual PH domains 1

Calculated MW:

374 aa, 43 kDa

Observed MW:

43 kDa

Purification Method:

Protein A purification

CloneNo.:

241370D6

Recommended Dilutions:

WB 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse brain tissue, HEK-293T cells, mouse small intestine tissue, rat kidney tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 84443-1-RR (ADAP1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



HEK-293T cells were subjected to SDS PAGE followed by western blot with 84443-1-RR (ADAP1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.