

For Research Use Only

# CD70 Recombinant antibody, PBS Only

Catalog Number: 84455-1-PBS



## Basic Information

Catalog Number:

84455-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC000725

GeneID (NCBI):

970

UNIPROT ID:

P32970

Full Name:

CD70 molecule

Calculated MW:

21 kDa

Observed MW:

28 kDa

Purification Method:

Protein A purification

CloneNo.:

241325D9

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



SKOV-3 cells were subjected to SDS PAGE followed by western blot with 84455-1-RR (CD70 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 84455-1-PBS in a different storage buffer formulation.