

For Research Use Only

Nectin-2/CD112 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84460-1-PBS



Basic Information

Catalog Number:

84460-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001042724.2

GeneID (NCBI):

5819

UNIPROT ID:

Q92692-2

Full Name:

poliovirus receptor-related 2
(herpesvirus entry mediator B)

Calculated MW:

51kDa

Purification Method:

Protein A purification

CloneNo.:

241725A10

Applications

Tested Applications:

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01329-1, Nectin-2/CD112 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84460-2-PBS. Detection antibody: 84460-1-PBS. Standard: Eg1670. Range: 0.781-50 ng/mL.



Cytometric bead array standard curve of MP01329-3, Nectin-2/CD112 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84460-3-PBS. Detection antibody: 84460-1-PBS. Standard: Eg1670. Range: 0.781-50 ng/mL.