

For Research Use Only

# Nectin-2/CD112 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84460-3-PBS



## Basic Information

Catalog Number:

84460-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_001042724.2

GeneID (NCBI):

5819

UNIPROT ID:

Q92692-2

Full Name:

poliovirus receptor-related 2  
(herpesvirus entry mediator B)

Calculated MW:

51kDa

Purification Method:

Protein A purification

CloneNo.:

241725C7

## Applications

Tested Applications:

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP01329-2, Nectin-2/CD112 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84460-3-PBS. Detection antibody: 84460-4-PBS. Standard: Eg1670. Range: 0.781-50 ng/mL



Cytometric bead array standard curve of MP01329-3, Nectin-2/CD112 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84460-3-PBS. Detection antibody: 84460-1-PBS. Standard: Eg1670. Range: 0.781-50 ng/mL