

For Research Use Only

CLEC4M Recombinant antibody, PBS Only (Detector)

Catalog Number: 84461-4-PBS



Basic Information

Catalog Number: 84461-4-PBS	GenBank Accession Number: BC038851	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 10332	CloneNo.: 241683E11
Source: Rabbit	UNIPROT ID: Q9H2X3	
Isotype: IgG	Full Name: C-type lectin domain family 4, member M	
	Calculated MW: 45 kDa	

Applications

Tested Applications:
Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01330-2, CLEC4M/DC-SIGNR Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84461-2-PBS. Detection antibody: 84461-4-PBS. Standard: RP00011. Range: 0.313-40 ng/mL.