

For Research Use Only

TLR2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84462-1-PBS



Basic Information

Catalog Number:

84462-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC033756

GeneID (NCBI):

7097

UNIPROT ID:

O60603

Full Name:

toll-like receptor 2

Calculated MW:

784 aa, 90 kDa

Purification Method:

Protein A purification

CloneNo.:

241680E6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01332-1, TLR2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84462-1-PBS. Detection antibody: 84462-4-PBS. Standard: RP00009. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP01332-2, TLR2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84462-1-PBS. Detection antibody: 84462-3-PBS. Standard: RP00009. Range: 0.781-100 ng/mL.