

For Research Use Only

# Integrin alpha-6 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84463-1-PBS



## Basic Information

Catalog Number: 84463-1-PBS	GenBank Accession Number: BC136456	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 3655	CloneNo.: 241703A12
Source: Rabbit	UNIPROT ID: P23229	
Isotype: IgG	Full Name: integrin, alpha 6	
	Calculated MW: 1073 aa, 119 kDa	

## Applications

Tested Applications:  
Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01334-1, Integrin alpha-6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84463-1-PBS. Detection antibody: 84463-2-PBS. Standard: RP00013. Range: 0.781-100 ng/mL