

For Research Use Only

LEPR Recombinant antibody, PBS Only (Detector)

Catalog Number: 84464-1-PBS



Basic Information

Catalog Number: 84464-1-PBS	GenBank Accession Number: NM_002303	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 3953	CloneNo.: 241692E8
Source: Rabbit	UNIPROT ID: P48357	
Isotype: IgG	Full Name: leptin receptor	
	Calculated MW: 132 kDa	

Applications

Tested Applications:
Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01335-1, Leptin R Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84464-3-PBS. Detection antibody: 84464-1-PBS. Standard: RP00010. Range: 1.56-100 ng/mL.