

For Research Use Only

# LEPR Recombinant antibody, PBS Only (Capture)

Catalog Number: 84464-3-PBS



## Basic Information

Catalog Number:

84464-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_002303

GeneID (NCBI):

3953

UNIPROT ID:

P48357

Full Name:

leptin receptor

Calculated MW:

132 kDa

Purification Method:

Protein A purification

CloneNo.:

241693C10

## Applications

Tested Applications:

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01335-1, Leptin R Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84464-3-PBS. Detection antibody: 84464-1-PBS. Standard: RP00010. Range: 1.56-100 ng/mL.



Cytometric bead array standard curve of MP01335-2, Leptin R Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84464-3-PBS. Detection antibody: 84464-2-PBS. Standard: RP00010. Range: 1.56-100 ng/mL.