

For Research Use Only

Integrin alpha-5 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84468-3-PBS



Basic Information

Catalog Number:

84468-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC008786

GeneID (NCBI):

3678

UNIPROT ID:

P08648

Full Name:

integrin, alpha 5 (fibronectin receptor,
alpha polypeptide)

Calculated MW:

115 kDa

Purification Method:

Protein A purification

CloneNo.:

241684G8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01350-1, Integrin alpha-5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84468-1-PBS. Detection antibody: 84468-3-PBS. Standard: RP00012. Range: 0.781-100 ng/mL