

For Research Use Only

Mouse Tie-2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84469-2-PBS



Basic Information

Catalog Number: 84469-2-PBS	GenBank Accession Number: NM_001290549.1	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 21687	CloneNo.: 241753E6
Source: Rabbit	UNIPROT ID: Q02858	
Isotype: IgG	Full Name: endothelial-specific receptor tyrosine kinase	
	Calculated MW: 126kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01351-1, MOUSE Tie-2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84469-2-PBS. Detection antibody: 84469-1-PBS. Standard: Eg1514. Range: 0.625-80 ng/mL.



Cytometric bead array standard curve of MP01351-3, MOUSE Tie-2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84469-2-PBS. Detection antibody: 84469-3-PBS. Standard: Eg1514. Range: 0.625-40 ng/mL.