

For Research Use Only

Integrin beta-6 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84475-2-PBS



Basic Information

Catalog Number: 84475-2-PBS	GenBank Accession Number: BC121178	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 3694	CloneNo.: 241877D8
Source: Rabbit	UNIPROT ID: P18564	
Isotype: IgG	Full Name: integrin, beta 6	
Immunogen Catalog Number: AG29002	Calculated MW: 788 aa, 86 kDa	

Applications

Tested Applications:
Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01331-1, Integrin beta-6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84475-2-PBS. Detection antibody: 84475-4-PBS. Standard: Ag29002. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP01331-2, Integrin beta-6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84475-2-PBS. Detection antibody: 84475-3-PBS. Standard: Ag29002. Range: 0.781-100 ng/mL