

For Research Use Only

ADAMTS2 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 84476-4-PBS



Basic Information

Catalog Number:

84476-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34957

GenBank Accession Number:

NM_014244

GeneID (NCBI):

9509

UNIPROT ID:

O95450

Full Name:

ADAM metalloproteinase with thrombospondin type 1 motif, 2

Calculated MW:

135 kDa

Purification Method:

Protein A purification

CloneNo.:

241860G5

Applications

Tested Applications:

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP01333-1, ADAMTS2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84476-2-PBS. Detection antibody: 84476-4-PBS. Standard: Ag34957. Range: 0.469-60 ng/mL.



Cytometric bead array standard curve of MP01333-2, ADAMTS2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84476-4-PBS. Detection antibody: 84476-3-PBS. Standard: Ag34957. Range: 0.469-60 ng/mL.



Cytometric bead array standard curve of MP01333-3, ADAMTS2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84476-4-PBS. Detection antibody: 84476-1-PBS. Standard: Ag34957. Range: 0.469-60 ng/mL.