

For Research Use Only

# CSF1R Recombinant antibody, PBS Only (Detector)

Catalog Number: 84477-1-PBS



## Basic Information

Catalog Number:

84477-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC047521

GeneID (NCBI):

1436

UNIPROT ID:

P07333

Full Name:

colony stimulating factor 1 receptor

Calculated MW:

972 aa, 108 kDa

Purification Method:

Protein A purification

CloneNo.:

241675B4

## Applications

Tested Applications:

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01340-1, CSF1R Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84477-3-PBS. Detection antibody: 84477-1-PBS. Standard: RP00001. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP01340-2, CSF1R Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84477-4-PBS. Detection antibody: 84477-1-PBS. Standard: RP00001. Range: 0.313-40 ng/mL