

For Research Use Only

CSF1R Recombinant antibody, PBS Only (Capture)

Catalog Number: 84477-4-PBS



Basic Information

Catalog Number:

84477-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC047521

GeneID (NCBI):

1436

UNIPROT ID:

P07333

Full Name:

colony stimulating factor 1 receptor

Calculated MW:

972 aa, 108 kDa

Purification Method:

Protein A purification

CloneNo.:

241675G3

Applications

Tested Applications:

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01340-2, CSF1R Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84477-4-PBS. Detection antibody: 84477-1-PBS. Standard: RP00001. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP01340-3, CSF1R Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84477-4-PBS. Detection antibody: 84477-2-PBS. Standard: CSF1R. Range: 0.313-40 ng/mL