

For Research Use Only

# GSR Recombinant antibody, PBS Only (Capture)

Catalog Number: 84480-2-PBS



## Basic Information

<b>Catalog Number:</b> 84480-2-PBS	<b>GenBank Accession Number:</b> BC069244	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 2936	<b>CloneNo.:</b> 241934D2
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> P00390	
<b>Isotype:</b> IgG	<b>Full Name:</b> glutathione reductase	
<b>Immunogen Catalog Number:</b> AG13080	<b>Calculated MW:</b> 522 aa, 56 kDa	

## Applications

**Tested Applications:**  
**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01316-1, GSR Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84480-2-PBS. Detection antibody: 84480-1-PBS. Standard: Ag13080. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP01316-2, GSR Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84480-2-PBS. Detection antibody: 84480-3-PBS. Standard: Ag13080. Range: 0.625-80 ng/mL