

For Research Use Only

# MCSF Recombinant antibody, PBS Only (Detector)

Catalog Number: 84481-1-PBS



## Basic Information

Catalog Number:

84481-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC021117

GeneID (NCBI):

1435

UNIPROT ID:

P09603

Full Name:

colony stimulating factor 1  
(macrophage)

Calculated MW:

60 kDa

Purification Method:

Protein A purification

CloneNo.:

241743A2

## Applications

Tested Applications:

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP01322-1, M-CSF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84481-3-PBS. Detection antibody: 84481-1-PBS. Standard: Eg1018. Range: 0.469-60 ng/mL