

For Research Use Only

# MCSF Recombinant antibody, PBS Only (Capture)

Catalog Number: 84481-3-PBS



## Basic Information

Catalog Number: 84481-3-PBS	GenBank Accession Number: BC021117	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1435	CloneNo.: 241743E6
Source: Rabbit	UNIPROT ID: P09603	
Isotype: IgG	Full Name: colony stimulating factor 1 (macrophage)	
	Calculated MW: 60 kDa	

## Applications

Tested Applications:  
Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01322-1, M-CSF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84481-3-PBS. Detection antibody: 84481-1-PBS. Standard: Eg1018. Range: 0.469-60 ng/mL



Cytometric bead array standard curve of MP01322-2, M-CSF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84481-3-PBS. Detection antibody: 84481-2-PBS. Standard: Eg1018. Range: 0.469-60 ng/mL