

For Research Use Only

HARBI1 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number:84482-2-PBS



Basic Information

Catalog Number: 84482-2-PBS	GenBank Accession Number: BC036925	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 283254	CloneNo.: 241915B11
Source: Rabbit	UNIPROT ID: Q96MB7	
Isotype: IgG	Full Name: harbinger transposase derived 1	
Immunogen Catalog Number: AG34137	Calculated MW: 349 aa, 39 kDa	

Applications

Tested Applications:
Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01317-1, HARBI 1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84482-2-PBS. Detection antibody: 84482-3-PBS. Standard: Ag34137. Range: 0.156-20 ng/mL



Cytometric bead array standard curve of MP01317-3, HARBI 1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84482-1-PBS. Detection antibody: 84482-2-PBS. Standard: Ag34137. Range: 0.156-20 ng/mL