

For Research Use Only

Mouse GM-CSF R alpha Recombinant antibody, PBS Only (Detector)

Catalog Number: 84483-3-PBS



Basic Information

Catalog Number:

84483-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_009970.2

GeneID (NCBI):

12982

UNIPROT ID:

Q00941

Full Name:

colony stimulating factor 2 receptor,
alpha, low-affinity (granulocyte-
macrophage)

Calculated MW:

42kDa

Purification Method:

Protein A purification

CloneNo.:

241889D7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01319-1, MOUSE CSF2RA/CD116 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84483-2-PBS. Detection antibody: 84483-3-PBS. Standard: Eg2016. Range: 0.469-60 ng/mL



Cytometric bead array standard curve of MP01319-3, MOUSE GM-CSF R alpha Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84483-1-PBS. Detection antibody: 84483-3-PBS. Standard: Eg2016. Range: 0.469-60 ng/mL