

For Research Use Only

CKAP4 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84485-3-PBS



Basic Information

Catalog Number: 84485-3-PBS	GenBank Accession Number: BC082972	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 10970	CloneNo.: 241871H5
Source: Rabbit	UNIPROT ID: Q07065	
Isotype: IgG	Full Name: cytoskeleton-associated protein 4	
Immunogen Catalog Number: AG10022	Calculated MW: 602 aa, 66 kDa	

Applications

Tested Applications:
Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01324-2, CKAP4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84485-1-PBS. Detection antibody: 84485-3-PBS. Standard: Ag10022. Range: 0.313-20 ng/mL