

For Research Use Only

GINS2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84486-1-PBS



Basic Information

Catalog Number:

84486-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9180

GenBank Accession Number:

BC010164

GeneID (NCBI):

51659

UNIPROT ID:

Q9Y248

Full Name:

GINS complex subunit 2 (Psf2 homolog)

Calculated MW:

185 aa, 21 kDa

Purification Method:

Protein A purification

CloneNo.:

241870B11

Applications

Tested Applications:

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01325-1, PSF2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84486-2-PBS. Detection antibody: 84486-1-PBS. Standard: Ag9180. Range: 0.313-20 ng/mL