

For Research Use Only

GIN52 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84486-3-PBS



Basic Information

Catalog Number: 84486-3-PBS	GenBank Accession Number: BC010164	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 51659	CloneNo.: 241870F6
Source: Rabbit	UNIPROT ID: Q9Y248	
Isotype: IgG	Full Name: GIN5 complex subunit 2 (Psf2 homolog)	
Immunogen Catalog Number: AG9180	Calculated MW: 185 aa, 21 kDa	

Applications

Tested Applications:
Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01325-2, PSF2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84486-2-PBS. Detection antibody: 84486-3-PBS. Standard: Ag9180. Range: 0.313-20 ng/mL