

For Research Use Only

Mouse EMMPRIN/CD147 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 84488-2-PBS



Basic Information

Catalog Number: 84488-2-PBS	GenBank Accession Number: NM_009768.2	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 12215	CloneNo.: 241852B9
Source: Rabbit	UNIPROT ID: P18572-1	
Isotype: IgG	Full Name: basigin	
	Calculated MW: 42kDa	

Applications

Tested Applications:
Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01320-1, MOUSE EMMPRIN/CD147 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84488-1-PBS. Detection antibody: 84488-4-PBS. Standard: Eg1530. Range: 0.781-50 ng/mL



Cytometric bead array standard curve of MP01320-2, MOUSE EMMPRIN/CD147 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84488-1-PBS. Detection antibody: 84488-2-PBS. Standard: Eg1530. Range: 0.781-50 ng/mL