

For Research Use Only

# Mouse EpCAM/TROP1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84490-1-PBS



## Basic Information

Catalog Number: 84490-1-PBS	GenBank Accession Number: NM_008532.2	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 17075	CloneNo.: 241851A12
Source: Rabbit	UNIPROT ID: Q99JW5	
Isotype: IgG	Full Name: epithelial cell adhesion molecule	
	Calculated MW: 35kDa	

## Applications

Tested Applications:  
Species Specificity:  
mouse

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01344-2, MOUSE EPCAM/CD326 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84490-2-PBS. Detection antibody: 84490-1-PBS. Standard: Eg1883. Range: 0.313-20 ng/mL



Cytometric bead array standard curve of MP01344-3, MOUSE EPCAM/CD326 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84490-3-PBS. Detection antibody: 84490-1-PBS. Standard: Eg1883. Range: 0.313-40 ng/mL