

For Research Use Only

Mouse EpCAM/TROP1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84490-2-PBS



Basic Information

Catalog Number:

84490-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_008532.2

GeneID (NCBI):

17075

UNIPROT ID:

Q99JW5

Full Name:

epithelial cell adhesion molecule

Calculated MW:

35kDa

Purification Method:

Protein A purification

CloneNo.:

241851E11

Applications

Tested Applications:

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01344-1, MOUSE EPCAM/CD326 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84490-2-PBS. Detection antibody: 84490-3-PBS. Standard: Eg1883. Range: 0.313-20 ng/mL



Cytometric bead array standard curve of MP01344-2, MOUSE EPCAM/CD326 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84490-2-PBS. Detection antibody: 84490-1-PBS. Standard: Eg1883. Range: 0.313-20 ng/mL