

For Research Use Only

PQBP1 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 84491-1-PBS



Basic Information

Catalog Number: 84491-1-PBS	GenBank Accession Number: BC012358	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 10084	CloneNo.: 241901D11
Source: Rabbit	UNIPROT ID: O60828	
Isotype: IgG	Full Name: polyglutamine binding protein 1	
Immunogen Catalog Number: AG9329	Calculated MW: 265 aa, 30 kDa	

Applications

Tested Applications:
Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01345-2, PQBP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84491-4-PBS. Detection antibody: 84491-1-PBS. Standard: Ag9329. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP01345-3, PQBP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84491-3-PBS. Detection antibody: 84491-1-PBS. Standard: Ag9329. Range: 0.313-40 ng/mL