

For Research Use Only

ACIII Recombinant antibody, PBS Only (Detector)

Catalog Number: 84492-1-PBS



Basic Information

Catalog Number: 84492-1-PBS	GenBank Accession Number: NM_004036	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 109	CloneNo.: 241906B5
Source: Rabbit	UNIPROT ID: O60266	
Isotype: IgG	Full Name: adenylate cyclase 3	
	Calculated MW: 129 kDa	

Applications

Tested Applications:
Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01346-2, ADCY3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84492-3-PBS. Detection antibody: 84492-1-PBS. Standard: Ag25859. Range: 1.563-100 ng/mL.